

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

A431 cells stimulated with EGF, Unconjugated NMB-X1003

Artikelname	A431 cells stimulated with EGF, Unconjugated
Artikelnummer	NMB-X1003
Hersteller Artikelnummer	X1003
Alternativnummer	NMB-X1003
Hersteller	NordicMubio
Kategorie	Proteine/Peptide
Applikation	WB
Konjugation	Unconjugated
Produktbeschreibung	A431 cells are a positive control for EGFR as they express high levels of EGFR. They also used for cell cycle studies and can be stimulated with EGF to stimulate EGFR and are therefore a very useful model to study EGFR downstream signalling cascades ...
Konzentration	See vial for concentration
Puffer	Provided as solution in Laemmli electrophoresis sample buffer. Samples need to be boiled for 3-5 minutes, before loading the gel.
Formulierung	Unconjugated
Formel	Provided as solution in Laemmli electrophoresis sample buffer. Samples need to be boiled for 3-5 minutes, before loading the gel.
Anwendungsbeschreibung	Suggested use: A volume of 10-25ul is suggested for immunoblotting.