

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Goat anti Bovine IgG (heavy and light chains), conjugated with Biotin, Clone: [Polyclonal], Monoclonal NMB-GAB/IGG(H+L)/BIO

Artikelname	Goat anti Bovine IgG (heavy and light chains), conjugated with Biotin, Clone: [Polyclonal], Monoclonal
Artikelnummer	NMB-GAB/IGG(H+L)/BIO
Hersteller Artikelnummer	GAB/IgG(H+L)/Bio
Alternativnummer	NMB-GAB/IGG(H+L)/BIO
Hersteller	NordicMubio
Wirt	Goat
Kategorie	Antikörper
Spezies Reaktivität	Bovine
Konjugation	Biotin
Format	IgG
Spezifität	IgG (H+L)
Minimale Kreuzreaktivität (MinX)	no cross-adsorbtion
Produktbeschreibung	The reactivity of the antiserum is directed to the Fc and Fab subunits of the IgG molecule. It includes a certain degree of reactivity with other immunoglobulins via the common Fab portion. The reactivity against Fab of IgG is less pronounced than in...
Klonalität	Monoclonal
Klon-Bezeichnung	[Polyclonal]

Puffer	Biotin-coupled purified hyperimmune goat IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added It is reconstituted by adding 1 ml sterile distilled water, spun down to remove insoluble particles, divided into sm
Quelle	Purified normal IgG isolated from pooled bovine serum. Freund's complete adjuvant is used in the first step of the immunization procedure.
Formulierung	Biotin
Formel	Biotin-coupled purified hyperimmune goat IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added
Antibody Type	Secondary Antibody
Anwendungsbeschreibung	ELISA, Immunocytochemistry, Immunohistochemistry (paraffin), Dot blot, Immunoblotting.