

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Goat anti Human lactoferrin, conjugated with Horseradish peroxidase, Clone: [Polyclonal], HRP, Monoclonal NMB-GAHU/LFR/PO

Artikelname	Goat anti Human lactoferrin, conjugated with Horseradish peroxidase, Clone: [Polyclonal], HRP, Monoclonal
Artikelnummer	NMB-GAHU/LFR/PO
Hersteller Artikelnummer	GAHu/Lfr/PO
Alternativnummer	NMB-GAHU/LFR/PO
Hersteller	NordicMubio
Wirt	Goat
Kategorie	Antikörper
Konjugation	HRP
Produktbeschreibung	In immunoelectrophoresis against human milk a single precipitin line is obtained. The antiserum does not react with any other protein component of human serum or plasma. In enzyme-immunocytochemical and immunohistochemical techniques for the detectio...
Klonalität	Monoclonal
Klon-Bezeichnung	[Polyclonal]
UniProt	P02788
Puffer	Horseradish peroxidase-coupled purified hyperimmune IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2) No preservative added, as it may interfere with the antibody activity.

Quelle	Exocrine organs produce various secretions, each with its characteristic function. Proteins found in secretions may be divided into two groups: those specific for the particular secretion, and plasma proteins independent of the type of exocrine cells. Lac
Formulierung	HRP
Formel	Horseradish peroxidase-coupled purified hyperimmune IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2) No preservative added, as it may interfere with the antibody activity.
Anwendungsbeschreibung	ELISA, Immunocytochemistry, Immunohistochemistry (paraffin), Dot blot, Immunoblotting. The lyophilized conjugate is shipped at ambient temperature and may be stored at +4C, prolonged storage at or below -20C. It is reconstituted by adding 1 ml sterile distilled water, spun down to remove insoluble particles, divided into small aliquots, frozen and stored at or below -20C. Prior to use, an aliquot is thawed slowly at ambient temperature, spun down again and used to prepare working dilutions by adding sterile phosphate buffered saline (PBS, pH 7.2). Repeated thawing and freezing should be avoided.