

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Goat anti Rat fibrinogen, conjugated with Horseradish peroxidase, Clone: [Polyclonal], HRP, Monoclonal NMB-GARA/FBG/PO

Artikelname	Goat anti Rat fibrinogen, conjugated with Horseradish peroxidase, Clone: [Polyclonal], HRP, Monoclonal
Artikelnummer	NMB-GARA/FBG/PO
Hersteller Artikelnummer	GARa/Fbg/PO
Alternativnummer	NMB-GARA/FBG/PO
Hersteller	NordicMubio
Wirt	Goat
Kategorie	Antikörper
Konjugation	HRP
Produktbeschreibung	The reactivity of the antiserum is restricted to fibrinogen. In immunoelectrophoresis and radial immunodiffusion (Ouchterlony), using various antiserum concentrations against normal rat plasma a single precipitin line is obtained which shows a reacti...
Klonalität	Monoclonal
Klon-Bezeichnung	[Polyclonal]
UniProt	P06399
Puffer	Horseradish peroxidase-coupled purified hyperimmune IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2) No preservative added, as it may interfere with the antibody activity.

Quelle	Fibrinogen (clotting factor I) is a heat labile beta glycoprotein present in plasma. It is the precursor of fibrin, which is the key protein constituting the network of the blood clot. Thrombin converts fibrinogen to fibrin by limited proteolysis. Fibrin
Formulierung	HRP
Formel	Horseradish peroxidase-coupled purified hyperimmune IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2) No preservative added, as it may interfere with the antibody activity.
Anwendungsbeschreibung	ELISA,Immunocytochemistry,Immunohistochemistry (paraffin),Dot blot,Immunoblotting. The lyophilized conjugate is shipped at ambient temperature and may be stored at +4C, prolonged storage at or below -20C. It is reconstituted by adding 1 ml sterile distilled water, spun down to remove insoluble particles, divided into small aliquots, frozen and stored at or below -20C. Prior to use, an aliquot is thawed slowly in the dark at ambient temperature, spun down again and used to prepare working dilutions by adding sterile phosphate buffered saline (PBS, pH 7.2). Repeated thawing and freezing should be avoided.