

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Rabbit anti Bacillus stearothermophilus Glycerokinase, conjugated with Biotin, Clone: [Polyclonal], Monoclonal NMB-NE069/BIO

Artikelname	Rabbit anti Bacillus stearothermophilus Glycerokinase, conjugated with Biotin, Clone: [Polyclonal], Monoclonal
Artikelnummer	NMB-NE069/BIO
Hersteller Artikelnummer	NE069/Bio
Alternativnummer	NMB-NE069/BIO
Hersteller	NordicMubio
Wirt	Rabbit
Kategorie	Antikörper
Konjugation	Biotin
Produktbeschreibung	The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and...
Klonalität	Monoclonal
Klon-Bezeichnung	[Polyclonal]
Puffer	Biotin-coupled hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added. It is reconstituted by adding 1.0 ml sterile dis

Quelle	Glycerokinase isolated and purified from <i>Bacillus stearothermophilus</i> . Freund's complete adjuvant is used in the first step of the immunization procedure.
Formulierung	Biotin
Formel	Biotin-coupled hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added.
Anwendungsbeschreibung	ELISA, Immunocytochemistry, Immunohistochemistry (paraffin), Dot blot, Immunoblotting.