

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Rabbit anti Human beta-2 microglobulin, Clone: [Polyclonal], Monoclonal NMB-RAHU/B2MG

Artikelname	Rabbit anti Human beta-2 microglobulin, Clone: [Polyclonal], Monoclonal
Artikelnummer	NMB-RAHU/B2MG
Hersteller Artikelnummer	RAHu/B2mg
Alternativnummer	NMB-RAHU/B2MG
Hersteller	NordicMubio
Wirt	Rabbit
Kategorie	Antikörper
Produktbeschreibung	The defined antibody reactivity is restricted to beta-2 microglobulin as tested at the level of sensitivity of immuno-precipitation techniques. A single precipitin line is obtained with foetal serum, colostrum and selected pathological plasma and uri...
Klonalität	Monoclonal
Klon-Bezeichnung	[Polyclonal]
UniProt	P61769
Puffer	Delipidated, heat inactivated, lyophilized, stable whole serum. No preservative added. Total protein and IgG concentrations in the antiserum are comparable to those of pooled normal rabbit serum. No foreign proteins added. Reconstitute the lyophilized ant

Quelle	Beta-2 microglobulin is a member of the Immunoglobulin superfamily. The amino acid sequence resembles that of immunoglobulin constant domains. It is isolated from urine of patients with impaired reabsorption in the kidney tubules. High levels of beta-2 mi
Formel	Delipidated, heat inactivated, lyophilized, stable whole serum. No preservative added. Total protein and IgG concentrations in the antiserum are comparable to those of pooled normal rabbit serum. No foreign proteins added.
Anwendungsbeschreibung	Precipitation assays. In immunoelectrophoresis use 2 µl serum, plasma or equivalent against 120 µl antiserum. In double radial immunodiffusion (Ouchterlony) use a rosette arrangement with 10 µl antiserum in a 3 mm diameter centre well and 2 µl serum samples (neat and serially diluted) in 2 mm diameter peripheral wells. In single radial immunodiffusion and electroimmunodiffusion use 0.5 to 1.0 percent antiserum in the agar gel.