

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Rabbit anti Mouse IgG1 IgG2a IgG2b IgG3 IgA IgM IgD IgE (heavy and light chains), Clone: [Polyclonal], Monoclonal NMB-RAM/IG(GAMDE)

Artikelname	Rabbit anti Mouse IgG1 IgG2a IgG2b IgG3 IgA IgM IgD IgE (heavy and light chains), Clone: [Polyclonal], Monoclonal
Artikelnummer	NMB-RAM/IG(GAMDE)
Hersteller Artikelnummer	RAM/Ig(GAMDE)
Alternativnummer	NMB-RAM/IG(GAMDE)
Hersteller	NordicMubio
Wirt	Rabbit
Kategorie	Antikörper
Spezies Reaktivität	Mouse
Konjugation	Unconjugated
Format	Antiserum
Spezifität	IgG+IgM+IgA+IgD+IgE (H+L)
Minimale Kreuzreaktivität (MinX)	no cross-adsorbtion
Produktbeschreibung	The reactivity of the antiserum is directed to the heavy chains of immunoglobulin isotypes IgG, IgA, IgM, IgD and IgE. The antiserum contains antibodies to common determinants and to the surface determinants of the common Fab as tested by immunoelect...
Klonalität	Monoclonal
Klon-Bezeichnung	[Polyclonal]

Puffer	Delipidated, heat inactivated, lyophilized stable whole serum. No preservative added, as it may interfere with the antibody activity. No foreign protein added. Total protein and IgG concentration in the antiserum are comparable to those of pooled rabbit s
Quelle	Purified polyclonal mouse IgG and pools of homogenous IgG, IgA, IgM, IgD and IgE. Freund's complete adjuvant is used in the first step of the immunization procedure.
Formel	Delipidated, heat inactivated, lyophilized stable whole serum. No preservative added, as it may interfere with the antibody activity. No foreign protein added. Total protein and IgG concentration in the antiserum are comparable to those of pooled rabbit
Antibody Type	Secondary Antibody
Anwendungsbeschreibung	Precipitation assays. In immunoelectrophoresis use 2 µl or equivalent against 120 µl antiserum. In double radial immunodiffusion use a rosette arrangement with 10 µl antiserum in a 3 mm diameter centre well and 2 µl serum samples (neat and diluted) in 2 mm diameter peripheral wells.