

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **Rabbit anti Photobacterium fischeri NADH-FMN Oxidoreductase, Clone: [Polyclonal], Monoclonal NMB-NE177/7S**

|                          |                                                                                                                                                                                                                                                                 |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Rabbit anti Photobacterium fischeri NADH-FMN Oxidoreductase, Clone: [Polyclonal], Monoclonal                                                                                                                                                                    |
| Artikelnummer            | NMB-NE177/7S                                                                                                                                                                                                                                                    |
| Hersteller Artikelnummer | NE177/7S                                                                                                                                                                                                                                                        |
| Alternativnummer         | NMB-NE177/7S                                                                                                                                                                                                                                                    |
| Hersteller               | NordicMubio                                                                                                                                                                                                                                                     |
| Wirt                     | Rabbit                                                                                                                                                                                                                                                          |
| Kategorie                | Antikörper                                                                                                                                                                                                                                                      |
| Produktbeschreibung      | The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and...   |
| Klonalität               | Monoclonal                                                                                                                                                                                                                                                      |
| Klon-Bezeichnung         | [Polyclonal]                                                                                                                                                                                                                                                    |
| Puffer                   | Purified hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added. It is reconstituted by adding 1.0 ml sterile distilled |
| Quelle                   | NADH-FMN oxidoreductase isolated and purified from Photobacterium fischeri. Freund's complete adjuvant is used in the first step of the immunization procedure.                                                                                                 |

|                        |                                                                                                                                                                                                          |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Formel                 | Purified hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added. |
| Anwendungsbeschreibung | Indirect immunofluorescence,ELISA,Dot blot,Immunoblotting.                                                                                                                                               |