

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **Rabbit IgG, conjugated with Horseradish peroxidase, HRP NMB-RIGG/PO**

|                          |                                                                                                                        |
|--------------------------|------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Rabbit IgG, conjugated with Horseradish peroxidase, HRP                                                                |
| Artikelnummer            | NMB-RIGG/PO                                                                                                            |
| Hersteller Artikelnummer | RIgG/PO                                                                                                                |
| Alternativnummer         | NMB-RIGG/PO                                                                                                            |
| Hersteller               | NordicMubio                                                                                                            |
| Kategorie                | Antikörper                                                                                                             |
| Konjugation              | HRP                                                                                                                    |
| Produktbeschreibung      | As a substitute for a non-relevant antibody, labelled with peroxidase or an alternative for the rabbit PAP complex.... |
| Puffer                   | Horseradish peroxidase coupled purified rabbit IgG in phosphate buffered saline (pH7.2). No preservative added         |
| Formulierung             | HRP                                                                                                                    |
| Formel                   | Horseradish peroxidase coupled purified rabbit IgG in phosphate buffered saline (pH7.2). No preservative added         |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Anwendungsbeschreibung | <p>ELISA,Immunocytochemistry,Immunohistochemistry (paraffin),Dot blot,Immunoblotting. The conjugate is shipped at ambient temperature and may be stored at +4C, prolonged storage at -20C. Reconstitute the lyophilized antiserum by adding 1 ml sterile distilled water. Spin down to remove insoluble particles, divide into small aliquots, and store frozen at or below -20C. Prior to use an aliquot is thawed at ambient temperature, spun down again, and used to prepare working dilutions by adding sterile PBS pH7.2). Repeated thawing and freezing should be avoided. If a slight precipitation occurs upon storage, this should be removed by centrifugation. It will not affect the performance of the conjugate. Diluted conjugate should be stored at +4C, not refrozen, and preferably used the same day.</p> |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|