

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **cAMP complete ELISA kit ENZ-ADI-900-163A**

|                          |                                                                                                                                                                                                                                                               |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | cAMP complete ELISA kit                                                                                                                                                                                                                                       |
| Artikelnummer            | ENZ-ADI-900-163A                                                                                                                                                                                                                                              |
| Hersteller Artikelnummer | ADI-900-163A                                                                                                                                                                                                                                                  |
| Alternativnummer         | ENZ-ADI-900-163A-96                                                                                                                                                                                                                                           |
| Hersteller               | Enzo Life Sciences                                                                                                                                                                                                                                            |
| Kategorie                | Kits/Assays                                                                                                                                                                                                                                                   |
| Applikation              | ELISA                                                                                                                                                                                                                                                         |
| Spezies Reaktivität      | All                                                                                                                                                                                                                                                           |
| Produktbeschreibung      | This is a colorimetric competitive immunoassay kit for the quantitative determination of extracellular cAMP diluted in buffer or intracellular cAMP in cells or tissues lysed in 0.1M HCl. 0.1M HCl is used to stop endogenous phosphodiesterase activity ... |
| Sensitivitaet            | Non-acetylated: 0.49 pmol/ml (assay buffer) / 1.18 pmol/ml (HCl) (range (non-acetylated) 0.78 - 200 pmol/ml)Acetylated: 0.027 pmol/ml (assay buffer) / 0.006 pmol/ml (HCl) (range (acetylated) 0.078 - 20 pmol/ml)                                            |
| Target-Kategorie         | cAMP                                                                                                                                                                                                                                                          |
| Anwendungsbeschreibung   | For the quantitative determination of cAMP in cell lysates, culture supernatants, serum, saliva, and tissue from any species.                                                                                                                                 |