

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Caldesmon, HMW (h-Caldesmon) (Smooth Muscle Marker) (rCALD1/820), 0.2mg/mL, Clone: [rCALD1/820], Mouse, Monoclonal BOT-BNUB1897-100

Artikelname	Caldesmon, HMW (h-Caldesmon) (Smooth Muscle Marker) (rCALD1/820), 0.2mg/mL, Clone: [rCALD1/820], Mouse, Monoclonal
Artikelnummer	BOT-BNUB1897-100
Hersteller Artikelnummer	BNUB1897-100
Alternativnummer	BOT-BNUB1897-100-100UL
Hersteller	Biotium
Wirt	Mouse
Kategorie	Antikörper
Applikation	IHC, WB
Spezies Reaktivität	Human
Immunogen	Recombinant full-length human CALD1 protein
Produktbeschreibung	Caldesmon HMW is the high molecular weight variant of Caldesmon. Two closely related variants of human caldesmon have been identified which are different in their electrophoretic mobility and cellular distribution. The h-caldesmon variant (120-150 kD...
Klonalität	Monoclonal
Konzentration	0.2 mg/mL
Klon-Bezeichnung	[rCALD1/820]
Molekulargewicht	150 kDa

UniProt	Q05682
Puffer	PBS, 0.05% BSA, 0.05% azide
Quelle	Animal
Anwendungsbeschreibung	Higher concentration may be required for direct detection using primary antibody conjugates than for indirect detection with secondary antibody Does not react with rat Immunohistology (formalin): 0.25-0.5 ug/mL for 30 minutes at RT Staining of formalin-fixed tissues requires boiling tissue sections in 1 mM EDTA pH 7.5-8.5 for 10-20 minutes followed by cooling at RT for 20 minutes Optimal dilution for a specific application should be determined by user