

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Neurofilament, phospho (NF-H) (Neuronal Marker)(NE14), 1mg/mL, Clone: [NE14], Mouse, Monoclonal BOT-BNUM1253-50

Artikelname	Neurofilament, phospho (NF-H) (Neuronal Marker)(NE14), 1mg/mL, Clone: [NE14], Mouse, Monoclonal
Artikelnummer	BOT-BNUM1253-50
Hersteller Artikelnummer	BNUM1253-50
Alternativnummer	BOT-BNUM1253-50-50UL
Hersteller	Biotium
Wirt	Mouse
Kategorie	Antikörper
Applikation	FC, IHC, WB
Spezies Reaktivität	Human, Mouse
Immunogen	Crude neurofilament preparation from porcine spinal cord
Produktbeschreibung	This MAb reacts with a 200 kDa protein, identified as heavy sub-unit of neurofilaments (NF-H). It reacts specifically with the phosphorylated KSP/KEP segment at the C-terminus of the heavy subunit (NF-H) of neurofilaments. After dephosphorylation of ...
Klonalität	Monoclonal
Konzentration	1 mg/mL
Klon-Bezeichnung	[NE14]
Molekulargewicht	200 kDa

UniProt	P12036
Puffer	PBS, no BSA, no azide
Quelle	Animal
Anwendungsbeschreibung	Higher concentration may be required for direct detection using primary antibody conjugates than for indirect detection with secondary antibody Immunohistochemistry (formalin-fixed): 0.25-0.5 ug/mL for 30 minutes at RT Western blot: 1-2 ug/mL Flow cytometry: 0.5-1 ug/million cells Staining of formalin-fixed tissues requires boiling tissue sections in 10 mM citrate buffer, pH 6.0, for 10-20 minutes followed by cooling at RT for 20 minutes Optimal dilution for a specific application should be determined by user