

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Recombinant Human BIRC5, E. coli EKL-PEH0132

|                          |                                                                                                                                                                                                                                                                                    |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Recombinant Human BIRC5, E. coli                                                                                                                                                                                                                                                   |
| Artikelnummer            | EKL-PEH0132                                                                                                                                                                                                                                                                        |
| Hersteller Artikelnummer | PEH0132                                                                                                                                                                                                                                                                            |
| Alternativnummer         | EKL-PEH0132-10UG,EKL-PEH0132-50UG,EKL-PEH0132-500UG                                                                                                                                                                                                                                |
| Hersteller               | EnkiLife                                                                                                                                                                                                                                                                           |
| Wirt                     | E. coli                                                                                                                                                                                                                                                                            |
| Kategorie                | Proteine/Peptide                                                                                                                                                                                                                                                                   |
| Spezies Reaktivität      | Human                                                                                                                                                                                                                                                                              |
| Produktbeschreibung      | Baculoviral IAP Repeat-Containing Protein 5 (BIRC5) belongs to the IAP family. BIRC5 exists as a homodimer in the apo state and as a monomer in the CPC-bound state. BIRC5 contains one BIR repeat and is expressed only in fetal kidney and liver, and to...                      |
| Molekulargewicht         | 16.4 KDa                                                                                                                                                                                                                                                                           |
| NCBI                     | <a href="#">15392</a>                                                                                                                                                                                                                                                              |
| Reinheit                 | Greater than 95% as determined by reducing SDS-PAGE                                                                                                                                                                                                                                |
| Anwendungsbeschreibung   | Always centrifuge tubes before opening.Do not mix by vortex or pipetting.It is not recommended to reconstitute to a concentration less than 100µg/ml.Dissolve the lyophilized protein in distilled water.Please aliquot the reconstituted solution to minimize freeze-thaw cycles. |