

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Recombinant Human SHH (C24II), E. coli EKL-PEH1552

|                          |                                                                                                                                                                                                                                                                                    |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Recombinant Human SHH (C24II), E. coli                                                                                                                                                                                                                                             |
| Artikelnummer            | EKL-PEH1552                                                                                                                                                                                                                                                                        |
| Hersteller Artikelnummer | PEH1552                                                                                                                                                                                                                                                                            |
| Alternativnummer         | EKL-PEH1552-10UG,EKL-PEH1552-50UG,EKL-PEH1552-500UG                                                                                                                                                                                                                                |
| Hersteller               | EnkiLife                                                                                                                                                                                                                                                                           |
| Wirt                     | E. coli                                                                                                                                                                                                                                                                            |
| Kategorie                | Proteine/Peptide                                                                                                                                                                                                                                                                   |
| Spezies Reaktivität      | Human                                                                                                                                                                                                                                                                              |
| Produktbeschreibung      | Sonic Hedgehog Homolog (SHH) belongs to a three-protein family called Hedgehog. The other two family members are Indian Hedgehog (IHH) and Desert Hedgehog (DHH). Hedgehog proteins are key signaling molecules in embryonic development. SHH is expressed...                      |
| Molekulargewicht         | 19.8 KDa                                                                                                                                                                                                                                                                           |
| NCBI                     | <a href="#">15465</a>                                                                                                                                                                                                                                                              |
| Reinheit                 | Greater than 95% as determined by reducing SDS-PAGE                                                                                                                                                                                                                                |
| Anwendungsbeschreibung   | Always centrifuge tubes before opening.Do not mix by vortex or pipetting.It is not recommended to reconstitute to a concentration less than 100µg/ml.Dissolve the lyophilized protein in distilled water.Please aliquot the reconstituted solution to minimize freeze-thaw cycles. |