

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Recombinant *S. cerevisiae* TIM16, *E. coli* EKL-PEV1650

|                          |                                                                                                                                                                                                                                                                                        |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Recombinant <i>S. cerevisiae</i> TIM16, <i>E. coli</i>                                                                                                                                                                                                                                 |
| Artikelnummer            | EKL-PEV1650                                                                                                                                                                                                                                                                            |
| Hersteller Artikelnummer | PEV1650                                                                                                                                                                                                                                                                                |
| Alternativnummer         | EKL-PEV1650-10UG,EKL-PEV1650-50UG,EKL-PEV1650-500UG                                                                                                                                                                                                                                    |
| Hersteller               | EnkiLife                                                                                                                                                                                                                                                                               |
| Wirt                     | <i>E. coli</i>                                                                                                                                                                                                                                                                         |
| Kategorie                | Proteine/Peptide                                                                                                                                                                                                                                                                       |
| Spezies Reaktivität      | Yeast                                                                                                                                                                                                                                                                                  |
| Produktbeschreibung      | Mitochondrial import inner membrane translocase subunit TIM16 (TIM16) is an essential component of the PAM complex. PAM complex is required for the translocation of transit peptide-containing proteins from the inner membrane into the mitochondrial mat...                         |
| Molekulargewicht         | 7.9 KDa                                                                                                                                                                                                                                                                                |
| NCBI                     | <a href="#">42949</a>                                                                                                                                                                                                                                                                  |
| Reinheit                 | Greater than 95% as determined by reducing SDS-PAGE                                                                                                                                                                                                                                    |
| Anwendungsbeschreibung   | Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. |