

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Recombinant Human TIMP-4 (C-6His) EKL-PHH1656

|                          |                                                                                                                                                                                                                                                                                    |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Recombinant Human TIMP-4 (C-6His)                                                                                                                                                                                                                                                  |
| Artikelnummer            | EKL-PHH1656                                                                                                                                                                                                                                                                        |
| Hersteller Artikelnummer | PHH1656                                                                                                                                                                                                                                                                            |
| Alternativnummer         | EKL-PHH1656-10UG,EKL-PHH1656-50UG,EKL-PHH1656-500UG                                                                                                                                                                                                                                |
| Hersteller               | EnkiLife                                                                                                                                                                                                                                                                           |
| Wirt                     | Human                                                                                                                                                                                                                                                                              |
| Kategorie                | Proteine/Peptide                                                                                                                                                                                                                                                                   |
| Spezies Reaktivität      | Human                                                                                                                                                                                                                                                                              |
| Produktbeschreibung      | Metalloproteinase inhibitor 4 is an enzyme that in humans is encoded by the TIMP4 gene, belongs to the protease inhibitor I35 (TIMP) family. The protein complexes with metalloproteinases (such as collagenases) and irreversibly inactivates them by bin...                      |
| Molekulargewicht         | 23.5 KDa                                                                                                                                                                                                                                                                           |
| NCBI                     | <a href="#">99727</a>                                                                                                                                                                                                                                                              |
| Reinheit                 | Greater than 95% as determined by reducing SDS-PAGE                                                                                                                                                                                                                                |
| Anwendungsbeschreibung   | Always centrifuge tubes before opening.Do not mix by vortex or pipetting.It is not recommended to reconstitute to a concentration less than 100µg/ml.Dissolve the lyophilized protein in distilled water.Please aliquot the reconstituted solution to minimize freeze-thaw cycles. |