

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Recombinant Mouse Siglec-2 (C-6His), Human EKL-PHM2219

|                          |                                                                                                                                                                                                                                                                                    |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Recombinant Mouse Siglec-2 (C-6His), Human                                                                                                                                                                                                                                         |
| Artikelnummer            | EKL-PHM2219                                                                                                                                                                                                                                                                        |
| Hersteller Artikelnummer | PHM2219                                                                                                                                                                                                                                                                            |
| Alternativnummer         | EKL-PHM2219-10UG,EKL-PHM2219-50UG,EKL-PHM2219-500UG                                                                                                                                                                                                                                |
| Hersteller               | EnkiLife                                                                                                                                                                                                                                                                           |
| Wirt                     | Human                                                                                                                                                                                                                                                                              |
| Kategorie                | Proteine/Peptide                                                                                                                                                                                                                                                                   |
| Spezies Reaktivität      | Mouse                                                                                                                                                                                                                                                                              |
| Produktbeschreibung      | Siglecs (sialic acid binding Ig-like lectins) are I-type (Ig-type) lectins belonging to the Ig superfamily. They are characterized by an N-terminal Ig-like V-type domain which mediates sialic acid binding, followed by varying numbers of Ig-like C2-ty...                      |
| Molekulargewicht         | 77.3 KDa                                                                                                                                                                                                                                                                           |
| NCBI                     | <a href="#">02562</a>                                                                                                                                                                                                                                                              |
| Reinheit                 | Greater than 95% as determined by reducing SDS-PAGE                                                                                                                                                                                                                                |
| Anwendungsbeschreibung   | Always centrifuge tubes before opening.Do not mix by vortex or pipetting.It is not recommended to reconstitute to a concentration less than 100µg/ml.Dissolve the lyophilized protein in distilled water.Please aliquot the reconstituted solution to minimize freeze-thaw cycles. |