

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Recombinant Rhesus Macaque ICOS (C-Fc), Human EKL-PHV2023

|                          |                                                                                                                                                                                                                                                                                    |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Recombinant Rhesus Macaque ICOS (C-Fc), Human                                                                                                                                                                                                                                      |
| Artikelnummer            | EKL-PHV2023                                                                                                                                                                                                                                                                        |
| Hersteller Artikelnummer | PHV2023                                                                                                                                                                                                                                                                            |
| Alternativnummer         | EKL-PHV2023-10UG,EKL-PHV2023-50UG,EKL-PHV2023-500UG                                                                                                                                                                                                                                |
| Hersteller               | EnkiLife                                                                                                                                                                                                                                                                           |
| Wirt                     | Human                                                                                                                                                                                                                                                                              |
| Kategorie                | Proteine/Peptide                                                                                                                                                                                                                                                                   |
| Spezies Reaktivität      | Monkey                                                                                                                                                                                                                                                                             |
| Produktbeschreibung      | Inducible T-cell costimulator, also known as activation-inducible lymphocyte immunomediatory molecule, CD278, AILIM, CVID1 and ICOS, belongs to the CD28 and CTLA4 cell surface receptor family.. ICOS contains one Ig-like V-type domain and exsits as a ...                      |
| Molekulargewicht         | 40.8 KDa                                                                                                                                                                                                                                                                           |
| NCBI                     | <a href="#">9</a>                                                                                                                                                                                                                                                                  |
| Reinheit                 | Greater than 95% as determined by reducing SDS-PAGE                                                                                                                                                                                                                                |
| Anwendungsbeschreibung   | Always centrifuge tubes before opening.Do not mix by vortex or pipetting.It is not recommended to reconstitute to a concentration less than 100µg/ml.Dissolve the lyophilized protein in distilled water.Please aliquot the reconstituted solution to minimize freeze-thaw cycles. |