

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### GENLISA Human Piezo-Type Mechanosensitive Ion Channel Component 2 (PIEZO2) ELISA KBS-KBH0883

|                          |                                                                                                                                                                                                                                                               |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | GENLISA Human Piezo-Type Mechanosensitive Ion Channel Component 2 (PIEZO2) ELISA                                                                                                                                                                              |
| Artikelnummer            | KBS-KBH0883                                                                                                                                                                                                                                                   |
| Hersteller Artikelnummer | KBH0883                                                                                                                                                                                                                                                       |
| Alternativnummer         | KBS-KBH0883-1                                                                                                                                                                                                                                                 |
| Hersteller               | Krishgen Biosystems                                                                                                                                                                                                                                           |
| Kategorie                | Kits/Assays                                                                                                                                                                                                                                                   |
| Spezies Reaktivität      | Human                                                                                                                                                                                                                                                         |
| Produktbeschreibung      | The GENLISA Human Piezo-Type Mechanosensitive Ion Channel Component 2 (PIEZO2) ELISA includes features like: Ready to use protocol with break-apart wells for ease of use Standardisation and High Reproducibility Lot to Lot Consistency Accuracy and Pre... |
| Detektionsbereich        | 0.156 - 10 ng/ml                                                                                                                                                                                                                                              |
| Sensitivitaet            | 0.05 ng/ml                                                                                                                                                                                                                                                    |
| Proben                   | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids. For tear samples, please connect with us to discuss further.                                                                                            |
| Anwendungsbeschreibung   | Colorimetric, Absorbance measured at 450nm                                                                                                                                                                                                                    |