

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### GENLISA Human AP-2 Complex Subunit Beta (AP2B1) ELISA KBS-KBH1234

|                          |                                                                                                                                                                                                                                                               |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | GENLISA Human AP-2 Complex Subunit Beta (AP2B1) ELISA                                                                                                                                                                                                         |
| Artikelnummer            | KBS-KBH1234                                                                                                                                                                                                                                                   |
| Hersteller Artikelnummer | KBH1234                                                                                                                                                                                                                                                       |
| Alternativnummer         | KBS-KBH1234-1                                                                                                                                                                                                                                                 |
| Hersteller               | Krishgen Biosystems                                                                                                                                                                                                                                           |
| Kategorie                | Kits/Assays                                                                                                                                                                                                                                                   |
| Spezies Reaktivität      | Human                                                                                                                                                                                                                                                         |
| Produktbeschreibung      | The GENLISA Human AP-2 Complex Subunit Beta (AP2B1) ELISA includes features like: Ready to use protocol with break-apart wells for ease of use Standardisation and High Reproducibility Lot to Lot Consistency Accuracy and Precision Validated against se... |
| Detektionsbereich        | 0.313 - 20 ng/ml                                                                                                                                                                                                                                              |
| Sensitivitaet            | 0.15 ng/ml                                                                                                                                                                                                                                                    |
| Proben                   | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids. For tear samples, please connect with us to discuss further.                                                                                            |
| Anwendungsbeschreibung   | Colorimetric, Absorbance measured at 450nm                                                                                                                                                                                                                    |