

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

MabSelect SuRe pcc protein A resin CRM-MA-331C

Artikelname	MabSelect SuRe pcc protein A resin
Artikelnummer	CRM-MA-331C
Hersteller Artikelnummer	Ma-331C
Alternativnummer	CRM-MA-331C-25
Hersteller	Creative BioMart
Kategorie	Proteine/Peptide
Anwendungsbeschreibung	<p>The MabSelect family, comprising MabSelect, MabSelect Xtra, MabSelect SuRe, and MabSelect SuRe LX chromatography media, was developed for high-capacity MAb capture at process scale. To enable the high mass transfer required for continuous chromatography, MabSelect SuRe pcc was further developed from MabSelect SuRe LX. The high dynamic binding capacity (DBC) for MAbs at short residence time contributes to increased productivity by reducing total time required for column load operations. Comparative studies with MabSelect SuRe LX show that higher capacities can be achieved at shorter residence times with MabSelect SuRe pcc.</p>