

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **Streptavidin Protein-ALP, Alkaline Phosphatase conjugated Streptavidin, AP, E. coli ABS-STN-NA117-1MG**

|                          |                                                                                                                                       |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Streptavidin Protein-ALP, Alkaline Phosphatase conjugated Streptavidin, AP, E. coli                                                   |
| Artikelnummer            | ABS-STN-NA117-1MG                                                                                                                     |
| Hersteller Artikelnummer | STN-NA117-1mg                                                                                                                         |
| Alternativnummer         | ABS-STN-NA117-1MG                                                                                                                     |
| Hersteller               | AcroBiosystems                                                                                                                        |
| Wirt                     | E. coli                                                                                                                               |
| Kategorie                | Proteine/Peptide                                                                                                                      |
| Spezies Reaktivität      | Bacteria                                                                                                                              |
| Konjugation              | AP                                                                                                                                    |
| Produktbeschreibung      | Streptavidin is expressed from E. coli cells and conjugated with alkaline phosphatase under optimal conditions. - Proteine/Peptide... |
| Molekulargewicht         | 13.8 kDa                                                                                                                              |
| Tag                      | Native                                                                                                                                |
| Puffer                   | 20 mM Tris, 3 M NaCl, pH7.5, 1% BSA, 0.03% Proclin300                                                                                 |
| Formulierung             | Liquid                                                                                                                                |
| Target-Kategorie         | Streptavidin                                                                                                                          |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Anwendungsbeschreibung | <p>Recommended for use in ALP Chemiluminescence System such as CLEIA (0.1 µg/mL) or MPCLIA, immunohistochemistry, and western blot applications. Avoid using biotin-containing solutions as diluents. It is recommended that the reagent be titrated for optimal performance for each application. NOTE: Do not use skim milk as a blocking agent in the assay with streptavidin, since skim milk contains free biotin which will cause high backgrounds.</p> |
|------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|