

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

ELISA Buffer Set (Phosphate system, 96T) **ABS-EBS-002-96TESTS**

Artikelname	ELISA Buffer Set (Phosphate system, 96T)
Artikelnummer	ABS-EBS-002-96TESTS
Hersteller Artikelnummer	EBS-002-96tests
Alternativnummer	ABS-EBS-002-96TESTS
Hersteller	AcroBiosystems
Kategorie	Kits/Assays
Produktbeschreibung	ELISA Buffer Set (Phosphate system, 96T) is expressed from . - Kits/Assays...
Formulierung	Solid
Anwendungsbeschreibung	Coating buffer: Use coating buffer to dilute the sample you will coat, add 100µL/well to high-bind plate. Washing buffer: Dilute the washing buffer to 1x with distilled water before experiment, add 300µL/well for washing the coated plate. Blocking buffer: Dilute the blocking buffer to 1x with distilled water before experiment, add 300µL/well for blocking the coated plate. The samples and detection antibodies can be prepared by 1x blocking buffer which should be 4 times dilution. Substrate solution: Add 100µL/well to plate. Stop solution: Add 50µL/well to plate.