

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

CAD Antibody / Caldesmon / CALD1 (HMW), Clone: [h-CALD], Mouse, Monoclonal NSJ-V2940SAF-100UG

Artikelname	CAD Antibody / Caldesmon / CALD1 (HMW), Clone: [h-CALD], Mouse, Monoclonal
Artikelnummer	NSJ-V2940SAF-100UG
Hersteller Artikelnummer	V2940SAF-100UG
Alternativnummer	NSJ-V2940SAF-100UG
Hersteller	NSJ Bioreagents
Wirt	Mouse
Kategorie	Antikörper
Applikation	FACS, IF, IHC-P, WB
Spezies Reaktivität	Human
Immunogen	Crude human uterus extract was used as the immunogen for the CAD antibody.
Produktbeschreibung	Caldesmon is an actin- and myosin-binding protein implicated in the regulation of actomyosin interactions in smooth muscle and nonmuscle cells (could act as a bridge between myosin and actin filaments). Stimulates actin binding of tropomyosin which i...
Klonalität	Monoclonal
Klon-Bezeichnung	[h-CALD]
UniProt	Q05682

Reinheit	Protein G affinity chromatography
Formulierung	1 mg/ml in 1X PBS, BSA free, sodium azide free
Antibody Type	Primary Antibody
Application Verdünnung	Flow cytometry: 1-2ug/10 ⁶ cells, Immunofluorescence: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT, Western blot: 1-2ug/ml
Anwendungsbeschreibung	Optimal dilution of the CAD antibody should be determined by the researcher. 1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 10-20 min followed by cooling at RT for 20 min. 2. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.