

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **CheKine(TM) Micro Glucose Oxidase Activity (GOD) Assay Kit ASC-KTB1310**

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | CheKine(TM) Micro Glucose Oxidase Activity (GOD) Assay Kit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Artikelnummer            | ASC-KTB1310                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Hersteller Artikelnummer | KTB1310                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Alternativnummer         | ASC-KTB1310-96, ASC-KTB1310-480                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Hersteller               | Abbkine Scientific                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Kategorie                | Kits/Assays                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Produktbeschreibung      | CheKine(TM) Glucose Oxidase Activity Assay Kit provides a simple and easy colorimetric assay for detecting Glucose Oxidase Activity in samples....                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Detektionsbereich        | 1-40µM                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Sensitivitaet            | 1µM                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Tag                      | GOD                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Anwendungsbeschreibung   | CheKine(TM) Micro Glucose Oxidase Activity Assay Kit provides a simple and easy colorimetric assay for detecting Glucose Oxidase Activity in biological samples such as in serum or plasma, cells, urine, other body fluid, food, growth medium. In this assay, GOD activity is determined by a coupled enzyme assay, in which GOD oxidizes D-glucose resulting in the production of hydrogen peroxide (H <sub>2</sub> O <sub>2</sub> ) that reacts with chromogen in an acid condition, produce a product that can be measured at OD 580 nm. Therefore, the glucose oxidase activity present in the sample is proportional to the signal obtained. |