

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Human TRAF2 and NCK-interacting protein kinase (TNIK) ELISA Kit ASC-KTE60200

Artikelname	Human TRAF2 and NCK-interacting protein kinase (TNIK) ELISA Kit
Artikelnummer	ASC-KTE60200
Hersteller Artikelnummer	KTE60200
Alternativnummer	ASC-KTE60200-48, ASC-KTE60200-96
Hersteller	Abbkine Scientific
Kategorie	Kits/Assays
Applikation	ELISA
Spezies Reaktivität	Human
Produktbeschreibung	This Human TRAF2 and NCK-interacting protein kinase (TNIK) ELISA Kit employs a two-site sandwich ELISA to quantitate TNIK....
Detektionsbereich	Please inquire
Sensitivitaet	Please inquire
Tag	TNIK
NCBI	23043
UniProt	Q9UKE5
Proben	Cell culture supernatantsSerumPlasmaOther biological fluids

Anwendungsbeschreibung	This Human TRAF2 and NCK-interacting protein kinase (TNIK) ELISA Kit employs a two-site sandwich ELISA to quantitate TNIK in samples. An antibody specific for TNIK has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any TNIK present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TNIK is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TNIK bound in the initial step. The color development is stopped and the intensity of the color is measured.
------------------------	--