

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Human Antigen peptide transporter 2 (TAP2) ELISA Kit ASC-KTE60472

Artikelname	Human Antigen peptide transporter 2 (TAP2) ELISA Kit
Artikelnummer	ASC-KTE60472
Hersteller Artikelnummer	KTE60472
Alternativnummer	ASC-KTE60472-48, ASC-KTE60472-96
Hersteller	Abbkine Scientific
Kategorie	Kits/Assays
Applikation	ELISA
Spezies Reaktivität	Human
Produktbeschreibung	This Human Antigen peptide transporter 2 (TAP2) ELISA Kit employs a two-site sandwich ELISA to quantitate TAP2....
Detektionsbereich	Please inquire
Sensitivitaet	Please inquire
Tag	TAP2
NCBI	6891
UniProt	Q03519
Proben	Cell culture supernatantsSerumPlasmaOther biological fluids

Anwendungsbeschreibung	This Human Antigen peptide transporter 2 (TAP2) ELISA Kit employs a two-site sandwich ELISA to quantitate TAP2 in samples. An antibody specific for TAP2 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any TAP2 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TAP2 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TAP2 bound in the initial step. The color development is stopped and the intensity of the color is measured.
------------------------	---