

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Human Mothers Against Decapentaplegic Homolog 2 (SMAD2) ELISA Kit ASC-KTE60575

Artikelname	Human Mothers Against Decapentaplegic Homolog 2 (SMAD2) ELISA Kit
Artikelnummer	ASC-KTE60575
Hersteller Artikelnummer	KTE60575
Alternativnummer	ASC-KTE60575-48, ASC-KTE60575-96
Hersteller	Abbkine Scientific
Kategorie	Kits/Assays
Applikation	ELISA
Spezies Reaktivität	Human
Produktbeschreibung	This Human Mothers Against Decapentaplegic Homolog 2 (SMAD2) ELISA Kit employs a two-site sandwich ELISA to quantitate SMAD2....
Detektionsbereich	Please inquire
Sensitivitaet	Please inquire
Tag	SMAD2
NCBI	4087
UniProt	Q15796
Proben	Cell culture supernatantsSerumPlasmaOther biological fluids

Anwendungsbeschreibung	This Human Mothers Against Decapentaplegic Homolog 2 (SMAD2) ELISA Kit employs a two-site sandwich ELISA to quantitate SMAD2 in samples. An antibody specific for SMAD2 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any SMAD2 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for SMAD2 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of SMAD2 bound in the initial step. The color development is stopped and the intensity of the color is measured.
------------------------	---