

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Human Resistin-like Molecule beta (RELM-beta) ELISA Kit ASC-KTE60803

|                          |                                                                                                                                  |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Human Resistin-like Molecule beta (RELM-beta) ELISA Kit                                                                          |
| Artikelnummer            | ASC-KTE60803                                                                                                                     |
| Hersteller Artikelnummer | KTE60803                                                                                                                         |
| Alternativnummer         | ASC-KTE60803-48, ASC-KTE60803-96                                                                                                 |
| Hersteller               | Abbkine Scientific                                                                                                               |
| Kategorie                | Kits/Assays                                                                                                                      |
| Applikation              | ELISA                                                                                                                            |
| Spezies Reaktivität      | Human                                                                                                                            |
| Produktbeschreibung      | Human Resistin-like Molecule beta (RELM-beta) ELISA Kit employs a two-site sandwich ELISA to quantitate RELM-beta in samples.... |
| Detektionsbereich        | 1.5 µmol/L-24 µmol/L                                                                                                             |
| Sensitivitaet            | 0.1 µmol/L                                                                                                                       |
| Tag                      | RELM-beta                                                                                                                        |
| Proben                   | Cell culture supernatantsSerumPlasmaOther biological fluids                                                                      |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Anwendungsbeschreibung | <p>Human Resistin-like Molecule beta (RELM-beta) ELISA Kit employs a two-site sandwich ELISA to quantitate RELM-beta in samples. An antibody specific for RELM-beta has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any RELM-beta present is bound by the immobilized antibody. After removing any unbound substances, HRP-Conjugate Human RELM-beta detection antibody is added to the wells. Following a wash to remove any unbound HRP reagent, a Chromogen solution is added to the wells and color develops in proportion to the amount of RELM-beta bound in the initial step. The color development is stopped and the intensity of the color is measured.</p> |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|