

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Human Phosphatidylinositol 3,4,5-trisphosphate-dependent Rac exchanger 2 protein (PREX2) ELISA Kit ASC-KTE61169

Artikelname	Human Phosphatidylinositol 3,4,5-trisphosphate-dependent Rac exchanger 2 protein (PREX2) ELISA Kit
Artikelnummer	ASC-KTE61169
Hersteller Artikelnummer	KTE61169
Alternativnummer	ASC-KTE61169-48, ASC-KTE61169-96
Hersteller	Abbkine Scientific
Kategorie	Kits/Assays
Applikation	ELISA
Spezies Reaktivität	Human
Produktbeschreibung	This Human Phosphatidylinositol 3,4,5-trisphosphate-dependent Rac exchanger 2 protein (PREX2) ELISA Kit employs a two-site sandwich ELISA to quantitate PREX2....
Detektionsbereich	Please inquire
Sensitivitaet	Please inquire
Tag	PREX2
NCBI	80243
UniProt	Q70Z35
Proben	Cell culture supernatantsSerumPlasmaOther biological fluids

Anwendungsbeschreibung	This Human Phosphatidylinositol 3,4,5-trisphosphate-dependent Rac exchanger 2 protein (PREX2) ELISA Kit employs a two-site sandwich ELISA to quantitate PREX2 in samples. An antibody specific for PREX2 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any PREX2 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for PREX2 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of PREX2 bound in the initial step. The color development is stopped and the intensity of the color is measured.
------------------------	--