

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Human Methemoglobin (MHB) ELISA Kit ASC-KTE61625

|                          |                                                                                                  |
|--------------------------|--------------------------------------------------------------------------------------------------|
| Artikelname              | Human Methemoglobin (MHB) ELISA Kit                                                              |
| Artikelnummer            | ASC-KTE61625                                                                                     |
| Hersteller Artikelnummer | KTE61625                                                                                         |
| Alternativnummer         | ASC-KTE61625-48, ASC-KTE61625-96                                                                 |
| Hersteller               | Abbkine Scientific                                                                               |
| Kategorie                | Kits/Assays                                                                                      |
| Applikation              | ELISA                                                                                            |
| Spezies Reaktivität      | Human                                                                                            |
| Produktbeschreibung      | This Human Methemoglobin (MHB) ELISA Kit employs a two-site sandwich ELISA to quantitate MHB.... |
| Detektionsbereich        | Please inquire                                                                                   |
| Sensitivitaet            | Please inquire                                                                                   |
| Tag                      | MHB                                                                                              |
| Proben                   | Cell culture supernatantsSerumPlasmaOther biological fluids                                      |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Anwendungsbeschreibung | <p>This Human Methemoglobin (MHB) ELISA Kit employs a two-site sandwich ELISA to quantitate MHB in samples. An antibody specific for MHB has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any MHB present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for MHB is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of MHB bound in the initial step. The color development is stopped and the intensity of the color is measured.</p> |
|------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|