

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Human Complement fragment 4b (C4b) ELISA Kit ASC-KTE62582

|                          |                                                                                                           |
|--------------------------|-----------------------------------------------------------------------------------------------------------|
| Artikelname              | Human Complement fragment 4b (C4b) ELISA Kit                                                              |
| Artikelnummer            | ASC-KTE62582                                                                                              |
| Hersteller Artikelnummer | KTE62582                                                                                                  |
| Alternativnummer         | ASC-KTE62582-48, ASC-KTE62582-96                                                                          |
| Hersteller               | Abbkine Scientific                                                                                        |
| Kategorie                | Kits/Assays                                                                                               |
| Applikation              | ELISA                                                                                                     |
| Spezies Reaktivität      | Human                                                                                                     |
| Produktbeschreibung      | This Human Complement fragment 4b (C4b) ELISA Kit employs a two-site sandwich ELISA to quantitate C4b.... |
| Detektionsbereich        | Please inquire                                                                                            |
| Sensitivitaet            | Please inquire                                                                                            |
| Tag                      | C4b                                                                                                       |
| Proben                   | Cell culture supernatantsSerumPlasmaOther biological fluids                                               |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Anwendungsbeschreibung | <p>This Human Complement fragment 4b (C4b) ELISA Kit employs a two-site sandwich ELISA to quantitate C4b in samples. An antibody specific for C4b has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any C4b present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for C4b is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of C4b bound in the initial step. The color development is stopped and the intensity of the color is measured.</p> |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|