

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Mouse Vasorin (VASN) ELISA Kit ASC-KTE70030

Artikelname	Mouse Vasorin (VASN) ELISA Kit
Artikelnummer	ASC-KTE70030
Hersteller Artikelnummer	KTE70030
Alternativnummer	ASC-KTE70030-48, ASC-KTE70030-96
Hersteller	Abbkine Scientific
Kategorie	Kits/Assays
Applikation	ELISA
Spezies Reaktivität	Mouse
Produktbeschreibung	This Mouse Vasorin (VASN) ELISA Kit employs a two-site sandwich ELISA to quantitate VASN....
Detektionsbereich	Please inquire
Sensitivitaet	Please inquire
Tag	VASN
NCBI	246154
UniProt	Q9CZT5
Proben	Cell culture supernatantsSerumPlasmaOther biological fluids

Anwendungsbeschreibung	This Mouse Vasorin (VASN) ELISA Kit employs a two-site sandwich ELISA to quantitate VASN in samples. An antibody specific for VASN has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any VASN present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for VASN is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of VASN bound in the initial step. The color development is stopped and the intensity of the color is measured.
------------------------	---