

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Porcine Heat Shock 70kDa Protein 1A (HSPA1A) ELISA Kit RDD-RD-HSPA1A-P

|                          |                                                                                                                                                                                                                                              |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Porcine Heat Shock 70kDa Protein 1A (HSPA1A) ELISA Kit                                                                                                                                                                                       |
| Artikelnummer            | RDD-RD-HSPA1A-P                                                                                                                                                                                                                              |
| Hersteller Artikelnummer | RD-HSPA1A-p                                                                                                                                                                                                                                  |
| Alternativnummer         | RDD-RD-HSPA1A-P-96,RDD-RD-HSPA1A-P-48                                                                                                                                                                                                        |
| Hersteller               | reddot Biotech                                                                                                                                                                                                                               |
| Kategorie                | Kits/Assays                                                                                                                                                                                                                                  |
| Applikation              | ELISA                                                                                                                                                                                                                                        |
| Spezies Reaktivität      | Porcine                                                                                                                                                                                                                                      |
| Produktbeschreibung      | This Porcine Heat Shock 70kDa Protein 1A (HSPA1A) ELISA Kit is a precoated, Sandwich ELISA assay designed to detect and quantify the target analyte in serum, plasma and other biological fluids . This product is for Research Use Only.... |
| Detektionsbereich        | 0.781-50ng/mL                                                                                                                                                                                                                                |
| Sensitivitaet            | 0.35ng/mL                                                                                                                                                                                                                                    |
| Proben                   | serum, plasma and other biological fluids                                                                                                                                                                                                    |
| Target-Kategorie         | Porcine Heat Shock 70kDa Protein 1A (HSPA1A)                                                                                                                                                                                                 |
| Anwendungsbeschreibung   | Assay Length: 4.5 hours. Experimental Method: Sandwich                                                                                                                                                                                       |