

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **ZAP 70 (TCR Zeta Associated Protein kinase 70kD), Rabbit USB-Z0300-01C**

|                          |                                                                                                                                                                                                                                                               |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | ZAP 70 (TCR Zeta Associated Protein kinase 70kD), Rabbit                                                                                                                                                                                                      |
| Artikelnummer            | USB-Z0300-01C                                                                                                                                                                                                                                                 |
| Hersteller Artikelnummer | Z0300-01C                                                                                                                                                                                                                                                     |
| Alternativnummer         | USB-Z0300-01C-100                                                                                                                                                                                                                                             |
| Hersteller               | US Biological                                                                                                                                                                                                                                                 |
| Wirt                     | Rabbit                                                                                                                                                                                                                                                        |
| Kategorie                | Antikörper                                                                                                                                                                                                                                                    |
| Applikation              | ELISA, IHC                                                                                                                                                                                                                                                    |
| Immunogen                | Synthesized non-phosphopeptide derived from human Zap-70 around the phosphorylation site of tyrosine 493 (S-Y-YP-T-A)                                                                                                                                         |
| Produktbeschreibung      | Applications: Suitable for use in ELISA and Immunohistochemistry. Other applications not tested. Recommended Dilution: Immunohistochemistry: 1:50-1:100 Optimal dilutions to be determined by the researcher. Storage and Stability: For long-term storage... |
| Reinheit                 | Purified by immunoaffinity chromatography.                                                                                                                                                                                                                    |
| Formulierung             | Supplied as a liquid in PBS (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM sodium chloride, 0.02% sodium azide and 50% glycerol.                                                                                                             |