

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Lactobacillus casei 60 kDa chaperonin (groL), partial BYT-ORB1785255

Artikelname	Recombinant Lactobacillus casei 60 kDa chaperonin (groL), partial
Artikelnummer	BYT-ORB1785255
Hersteller Artikelnummer	orb1785255
Alternativnummer	BYT-ORB1785255-1,BYT-ORB1785255-100,BYT-ORB1785255-20
Hersteller	Biorbyt
Kategorie	Proteine/Peptide
Produktbeschreibung	This Recombinant Lactobacillus casei 60 kDa chaperonin (groL), partial spans the amino acid sequence from region 182-374aa. Purity: Greater than 90% as determined by SDS-PAGE....
Molekulargewicht	27.9 kDa
UniProt	B3W9W7
Puffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Quelle	Lacticaseibacillus casei (strain BL23) (Lactobacillus casei)
Reinheit	Greater than 90% as determined by SDS-PAGE.
Formulierung	Liquid or Lyophilized powder

Sequenz	DTELSVVEGMQFDRGYLSQYMVTDNDKMEADLDDPYILITDKKISNIQDILPLL QEIVQQGKALLIIADDVAGEALPTLVLNKIRGTFNVVAVKAPGFGDRRKAQLED IATLTGGTVISSDLGLDLKDTKLEQLGRAGKVTVTKDNTTIVDGAGSKDAIAER VNIKKQIDDTTSDFDREKLQERLAKLAGGV
Anwendungsbeschreibung	Biological Origin: Lacticaseibacillus casei (strain BL23) (Lactobacillus casei). Application Notes: We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference