

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Recombinant *Pseudomonas aeruginosa* Major outer membrane lipoprotein (oprI) BYT-ORB1785453

|                          |  |
|--------------------------|--|
| Artikelname              | Recombinant <i>Pseudomonas aeruginosa</i> Major outer membrane lipoprotein (oprI)  |
| Artikelnummer            | BYT-ORB1785453   |
| Hersteller Artikelnummer | orb1785453   |
| Alternativnummer         | BYT-ORB1785453-1,BYT-ORB1785453-100,BYT-ORB1785453-20  |
| Hersteller               | Biorbyt  |
| Kategorie                | Proteine/Peptide   |
| Produktbeschreibung      | This Recombinant <i>Pseudomonas aeruginosa</i> Major outer membrane lipoprotein (oprI) spans the amino acid sequence from region 1-83aa. Purity: Greater than 90% as determined by SDS-PAGE....                                    |
| Molekulargewicht         | 40.3 kDa   |
| UniProt                  | <a href="#">P11221</a>   |
| Puffer                   | If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0. |
| Quelle                   | <i>Pseudomonas aeruginosa</i> (strain ATCC 15692 / DSM 22644 / CIP 104116 / JCM 14847 / LMG 12228 / 1C / PRS 101 / PAO1)   |
| Reinheit                 | Greater than 90% as determined by SDS-PAGE.  |
| Formulierung             | Liquid or Lyophilized powder   |

|                        |   |
|------------------------|---|
| Sequenz                | MNNVLKFSALALAAVLATGCSSHSKETEARLTATEDAAARAQARADEAYRKA<br>DEALGAAQKAQQTADEANERALRMLEKASRK   |
| Anwendungsbeschreibung | Biological Origin: Pseudomonas aeruginosa (strain ATCC 15692 / DSM 22644 / CIP 104116 / JCM 14847 / LMG 12228 / 1C / PRS 101 / PAO1). Application Notes: We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference |