

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Sheep anti Guinea Pig Serum proteins, Clone: [Polyclonal], Monoclonal NMB-SHAGP/TSP

|                          |                                                                                                                                                                                                                                                                |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Sheep anti Guinea Pig Serum proteins, Clone: [Polyclonal], Monoclonal                                                                                                                                                                                          |
| Artikelnummer            | NMB-SHAGP/TSP                                                                                                                                                                                                                                                  |
| Hersteller Artikelnummer | ShAGp/TSP                                                                                                                                                                                                                                                      |
| Alternativnummer         | NMB-SHAGP/TSP                                                                                                                                                                                                                                                  |
| Hersteller               | NordicMubio                                                                                                                                                                                                                                                    |
| Wirt                     | Sheep                                                                                                                                                                                                                                                          |
| Kategorie                | Antikörper                                                                                                                                                                                                                                                     |
| Produktbeschreibung      | In immunoelectrophoresis against pooled serum and enriched serum proteins fractions precipitation can be observed of not less than 18 individual proteins components. In precipitating techniques as immunoelectrophoresis and radial immunodiffusion (Ouc...  |
| Klonalität               | Monoclonal                                                                                                                                                                                                                                                     |
| Klon-Bezeichnung         | [Polyclonal]                                                                                                                                                                                                                                                   |
| Puffer                   | Delipidated, heat inactivated, lyophilized, stable whole antiserum. No preservative added. Total protein and IgG concentrations in the antiserum are comparable to those of pooled normal sheep serum. No foreign proteins added. Reconstitute the lyophilized |
| Quelle                   | Pooled whole guinea pig serum and partly purified serum fractions. Freund's complete adjuvant is used in the first step of the immunization procedure.                                                                                                         |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Formel                 | Delipidated, heat inactivated, lyophilized, stable whole antiserum. No preservative added. Total protein and IgG concentrations in the antiserum are comparable to those of pooled normal sheep serum. No foreign proteins added.                                                                                                                                                                                                                                                                       |
| Anwendungsbeschreibung | <p>Precipitation assays. In immunoelectrophoresis use 2 µl serum or equivalent against 120 µl antiserum. In double radial immunodiffusion (Ouchterlony) use a rosette arrangement with 10 µl antiserum in 3 mm diameter center well and 2 µl serum samples (neat and serially diluted in 2 mm diameter peripheral wells. Different bleedings of the immunized animals are pooled to obtain a broad spectrum balanced against the varying concentrations of the individual serum protein components.</p> |