

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Swine anti Goat IgM (Fc specific), Clone: [Polyclonal], Porcine, Monoclonal NMB-SWAG/IGM(FC)

|                                  |                                                                                                                                                                                                                                                               |
|----------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname                      | Swine anti Goat IgM (Fc specific), Clone: [Polyclonal], Porcine, Monoclonal                                                                                                                                                                                   |
| Artikelnummer                    | NMB-SWAG/IGM(FC)                                                                                                                                                                                                                                              |
| Hersteller Artikelnummer         | SwAG/IgM(Fc)                                                                                                                                                                                                                                                  |
| Alternativnummer                 | NMB-SWAG/IGM(FC)                                                                                                                                                                                                                                              |
| Hersteller                       | NordicMubio                                                                                                                                                                                                                                                   |
| Wirt                             | Porcine                                                                                                                                                                                                                                                       |
| Kategorie                        | Antikörper                                                                                                                                                                                                                                                    |
| Spezies Reaktivität              | Goat                                                                                                                                                                                                                                                          |
| Konjugation                      | Unconjugated                                                                                                                                                                                                                                                  |
| Format                           | Antiserum                                                                                                                                                                                                                                                     |
| Spezifität                       | IgM (μ)                                                                                                                                                                                                                                                       |
| Minimale Kreuzreaktivität (MinX) | no cross-adsorbtion                                                                                                                                                                                                                                           |
| Produktbeschreibung              | The reactivity of the antiserum is restricted to the Fc part of the IgM molecule. In immunoelectrophoresis and radial immunodiffusion, using various antiserum concentrations against normal goat serum a single precipitin line is obtained which shows a... |
| Klonalität                       | Monoclonal                                                                                                                                                                                                                                                    |
| Klon-Bezeichnung                 | [Polyclonal]                                                                                                                                                                                                                                                  |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                |
|------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Puffer                 | Delipidated, heat inactivated, lyophilized, stable whole antiserum. No preservative added Total protein and IgG concentrations in the antiserum are comparable to those of pooled normal swine serum. No foreign proteins added. Reconstitute the lyophilized a                                                                                                                                                                                |
| Quelle                 | Highly purified normal IgM isolated from pooled goat serum. Freund's complete adjuvant is used in the first step of the immunization procedure.                                                                                                                                                                                                                                                                                                |
| Formel                 | Delipidated, heat inactivated, lyophilized, stable whole antiserum. No preservative added Total protein and IgG concentrations in the antiserum are comparable to those of pooled normal swine serum. No foreign proteins added.                                                                                                                                                                                                               |
| Antibody Type          | Secondary Antibody                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Anwendungsbeschreibung | Precipitation assays. In immunoelectrophoresis use 2 µl serum, plasma or equivalent against 120 µl antiserum. In double radial immunodiffusion (Ouchterlony) use a rosette arrangement with 10 µl antiserum in a 3 mm diameter centre well and 2 µl serum samples (neat and serially diluted) in 2 mm diameter peripheral wells. In single radial immunodiffusion and electroimmunodiffusion use 0.5 to 1.0 percent antiserum in the agar gel. |