

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Swine anti Mouse IgG (heavy and light chains), Clone: [Polyclonal], Porcine, Monoclonal NMB-SWAM/IGG(H+L)

Artikelname	Swine anti Mouse IgG (heavy and light chains), Clone: [Polyclonal], Porcine, Monoclonal
Artikelnummer	NMB-SWAM/IGG(H+L)
Hersteller Artikelnummer	SwAM/IgG(H+L)
Alternativnummer	NMB-SWAM/IGG(H+L)
Hersteller	NordicMubio
Wirt	Porcine
Kategorie	Antikörper
Applikation	IEP, IP
Spezies Reaktivität	Mouse
Konjugation	Unconjugated
Format	Antiserum
Spezifität	IgG (H+L)
Minimale Kreuzreaktivität (MinX)	no cross-adsorbtion
Produktbeschreibung	The reactivity of the antiserum is directed to the Fc and Fab subunits of the IgG molecule. It includes a certain degree of reactivity with other immunoglobulins via the common Fab portion. It does not react with any non-Ig protein in mouse serum, as...
Klonalität	Monoclonal

Klon-Bezeichnung	[Polyclonal]
Puffer	Neat liquid antiserum containing 0.05% sodium azide as a preservative.
Quelle	Purified normal IgG isolated from pooled mouse serum. Freund's complete adjuvant is used in the first step of the immunization procedure.
Formel	Neat liquid antiserum containing 0.05% sodium azide as a preservative.
Antibody Type	Secondary Antibody
Anwendungsbeschreibung	Precipitation assays. In immunoelectrophoresis use 2 µl or equivalent against 120 µl antiserum. In double radial immunodiffusion (Ouchterlony) use a rosette arrangement with 10 µl antiserum in a 3 mm diameter centre well and 2 µl serum samples (neat and diluted) in 2 mm diameter peripheral wells.