

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### ELISpot Pro: Human IFN-gamma (HRP) MAB-3420-2HPW-10

|                          |                                                                                                                                                                                                                                                               |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | ELISpot Pro: Human IFN-gamma (HRP)                                                                                                                                                                                                                            |
| Artikelnummer            | MAB-3420-2HPW-10                                                                                                                                                                                                                                              |
| Hersteller Artikelnummer | 3420-2HPW-10                                                                                                                                                                                                                                                  |
| Alternativnummer         | MAB-3420-2HPW-10                                                                                                                                                                                                                                              |
| Hersteller               | MABTECH                                                                                                                                                                                                                                                       |
| Kategorie                | Kits/Assays                                                                                                                                                                                                                                                   |
| Applikation              | ELISpot                                                                                                                                                                                                                                                       |
| Spezies Reaktivität      | Human                                                                                                                                                                                                                                                         |
| Konjugation              | HRP                                                                                                                                                                                                                                                           |
| Produktbeschreibung      | ELISpot Pro kit for enumeration of cells secreting IFN-gamma, The ELISpot Pro kit is Mabtechs most innovative product featuring pre-coated plates and one-step detection. The kit enables enumeration of cells secreting human IFN-gamma and contains ELIS... |
| Detektionsbereich        | One-step HRP                                                                                                                                                                                                                                                  |
| NCBI                     | <a href="#">3458</a>                                                                                                                                                                                                                                          |
| Pubmed                   | <a href="https://pubmed.ncbi.nlm.nih.gov/11792386">https://pubmed.ncbi.nlm.nih.gov/11792386</a> target=_blank rel=nofollow>11792386                                                                                                                           |
| Target-Kategorie         | IFN-gamma                                                                                                                                                                                                                                                     |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Anwendungsbeschreibung | <p>Complementary Products:<br/> <a href="https://www.mabtech.com/products/peppool-cefras-global-cd8-human-3637-1-0">https://www.mabtech.com/products/peppool-cefras-global-cd8-human-3637-1-0</a>, <a href="https://www.mabtech.com/products/peppool-cef-cd8-human-3616-1-0">https://www.mabtech.com/products/peppool-cef-cd8-human-3616-1-0</a>, <a href="https://www.mabtech.com/products/peppool-cefta-cd4-human-3617-1-0">https://www.mabtech.com/products/peppool-cefta-cd4-human-3617-1-0</a>,<br/> <a href="https://www.mabtech.com/products/peppool-cefras-global-cd4-human-3627-1-0">https://www.mabtech.com/products/peppool-cefras-global-cd4-human-3627-1-0</a>, <a href="https://www.mabtech.com/products/peppool-cefras-global-cd4-and-cd8-human-3638-1-0">https://www.mabtech.com/products/peppool-cefras-global-cd4-and-cd8-human-3638-1-0</a>,<br/> <a href="https://www.mabtech.com/products/peppool-cefras-global-cd8-human-3637-1-0">https://www.mabtech.com/products/peppool-cefras-global-cd8-human-3637-1-0</a>, <a href="https://www.mabtech.com/products/peppool-cefras-global-cd4-and-cd8-human-3638-1-0">https://www.mabtech.com/products/peppool-cefras-global-cd4-and-cd8-human-3638-1-0</a>. Contents: 7-B6-1, HRP,TMB substrate for ELISpot,Pre-coated ELISpot plate, white MSIP (mAb 1-D1K),Anti-CD3 mAb (CD3-2)</p> |
|------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|