

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

ELISpot Plus: Human IFN-gamma (ALP) MAB-3420-4AST-2

Artikelname	ELISpot Plus: Human IFN-gamma (ALP)
Artikelnummer	MAB-3420-4AST-2
Hersteller Artikelnummer	3420-4AST-2
Alternativnummer	MAB-3420-4AST-2
Hersteller	MABTECH
Kategorie	Kits/Assays
Applikation	ELISpot
Spezies Reaktivität	Human
Konjugation	ALP
Produktbeschreibung	ELISpot Plus for enumeration of cells secreting IFN-gamma, This kit is ideal for users who want a convenient and sensitive assay. The assay is designed for the enumeration of cells secreting human IFN-gamma. The kit includes 96-well ELISpot strip pla...
Detektionsbereich	Streptavidin-ALP
NCBI	3458
Pubmed	https://pubmed.ncbi.nlm.nih.gov/11792386 target=_blank rel=nofollow>11792386
Target-Kategorie	IFN-gamma

Anwendungsbeschreibung	Complementary Products: https://www.mabtech.com/products/peppool-cefras-global-cd8-human-3637-1-0, https://www.mabtech.com/products/peppool-cef-cd8-human-3616-1-0, https://www.mabtech.com/products/peppool-cefta-cd4-human-3617-1-0, https://www.mabtech.com/products/peppool-cefras-global-cd4-human-3627-1-0, https://www.mabtech.com/products/peppool-cefras-global-cd4-and-cd8-human-3638-1-0, https://www.mabtech.com/products/peppool-cefras-global-cd8-human-3637-1-0, https://www.mabtech.com/products/peppool-cefras-global-cd4-and-cd8-human-3638-1-0. Contents: 7-B6-1, biotin,Streptavidin-ALP enzyme conjugate for ELISpot,BCIP/NBT-plus substrate for ELISpot,Anti-CD3 mAb (CD3-2),Pre-coated ELISpot plate, strip (mAb 1-D1K)
------------------------	--