

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Goat IgG anti-Mouse IgG (H+L)-unconj., MinX none DNA-SEC-183158

Artikelname	Goat IgG anti-Mouse IgG (H+L)-unconj., MinX none
Artikelnummer	DNA-SEC-183158
Hersteller Artikelnummer	SEC-183158
Alternativnummer	DNA-SEC-183158
Hersteller	dianova
Wirt	Goat
Kategorie	Antikörper
Applikation	ELISA,FLISA,FACS,IHC,IF,IP,WB
Spezies Reaktivität	Mouse
Immunogen	Anti-Mouse IgG whole molecule was produced by repeated immunization with Mouse IgG whole molecule in goat.
Konjugation	Unconjugated
Format	IgG
Spezifität	IgG (H+L)
Minimale Kreuzreaktivität (MinX)	no cross-adsorbtion
Produktbeschreibung	Anti-Mouse IgG Antibody generated in goat detects reactivity to Mouse IgG. Secreted as part of the adaptive immune response by plasma B cells, immunoglobulin G constitutes 75% of serum immunoglobulins. Immunoglobulin G binds to viruses, bacteria, as ...

Klonalität	Polyclonal
Konzentration	2.19 mg/mL
Isotyp	Ig
Puffer	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Reinheit	Secondary Antibody Anti-Mouse IgG (H&L) was prepared from monospecific antiserum by immunoaffinity chromatography using Mouse IgG coupled to agarose. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Goat Serum, Mouse IgG and Mouse Serum.
Formulierung	Liquid (sterile filtered)
Formel	20 mM K3PO4,150 mM NaCl,pH 7,2,sterile filtered,0,01% NaN3
Target-Kategorie	Mouse
Antibody Type	Secondary Antibody
Application Verdünnung	ELISA Dilution: 1:25,000, FLISA Dilution: User Optimized, Flow Cytometry Dilution: User Optimized, Immunohistochemistry Dilution: 1:1,000 - 1:5,000, IF Microscopy Dilution: User Optimized, Immunoprecipitation Dilution: User Optimized, Western Blot Diluti
Anwendungsbeschreibung	Anti-Mouse IgG affinity purified secondary antibody is generated in goat detects specifically Mouse IgG whole molecule. This anti-Mouse IgG secondary antibody is ideal for investigators who routinely perform ELISA, Sandwich ELISA, titration assays, western-blot, immunoprecipitation and more generally immunoassays. Specific conditions for reactivity and signal detection should be optimized by the end user.