

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Rabbit IgG anti-Sheep IgG (H+L)-HRPO, MinX Hu DNA-SEC-183552

Artikelname	Rabbit IgG anti-Sheep IgG (H+L)-HRPO, MinX Hu
Artikelnummer	DNA-SEC-183552
Hersteller Artikelnummer	SEC-183552
Alternativnummer	DNA-SEC-183552
Hersteller	dianova
Wirt	Rabbit
Kategorie	Antikörper
Applikation	ELISA,IHC,WB
Spezies Reaktivität	Sheep
Immunogen	Sheep IgG whole molecule
Konjugation	HRPO
Format	IgG
Spezifität	IgG (H+L)
Minimale Kreuzreaktivität (MinX)	Human
Produktbeschreibung	Anti-Sheep IgG generated in rabbit detects sheep Immunoglobulin G. It is a protein complex composed of four peptide chains - two identical heavy chains and two identical light chains arranged in a Y-shape typical of antibody monomers. Each IgG has tw...
Klonalität	Polyclonal

Konzentration	1.0 mg/mL
Isotyp	Ig
Puffer	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Reinheit	Anti-Sheep IgG HRP Conjugated Antibody Pre-Adsorbed was prepared from monospecific antiserum by immunoaffinity chromatography using Sheep IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Peroxidase, anti-Rabbit Serum, Sheep IgG and Sheep Serum. No reaction was observed against Human Serum Proteins.
Formulierung	Lyophilized
Formel	20 mM K3PO4,150 mM NaCl,pH 7,2,lyophilisate,0,01% Gentamicin
Target-Kategorie	Sheep
Antibody Type	Secondary Antibody
Application Verdünnung	ELISA Dilution: 1:30,000 - 1:70,000, Immunohistochemistry Dilution: 1:250 - 1:1,500, Western Blot Dilution: 1:1,000 - 1:10,000
Anwendungsbeschreibung	Anti-Sheep IgG Antibody Pre-Adsorbed conjugated to horseradish peroxidase is available in a variety of formats. HRP Anti-Sheep IgG Secondary Antibody conjugate is suitable for western blot or dot blot, ELISA and immunohistochemistry, immunoperoxidase electron microscopy and as well as other enzymatic assays requiring lot-to-lot consistency.