

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### 4-component ICP-MS Tuning Solution, Mix 2 at 10 µg/mL 100 mL Contains: HPS-ICP-MS-TS-2-100

|                            |                                                                                                                                                                                                                                                 |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | 4-component ICP-MS Tuning Solution, Mix 2 at 10 µg/mL 100 mL Contains:                                                                                                                                                                          |
| Biozol Catalog Number      | HPS-ICP-MS-TS-2-100                                                                                                                                                                                                                             |
| Supplier Catalog Number    | ICP-MS-TS-2-100                                                                                                                                                                                                                                 |
| Alternative Catalog Number | HPS-ICP-MS-TS-2-100-100                                                                                                                                                                                                                         |
| Manufacturer               | High Purity Standards                                                                                                                                                                                                                           |
| Category                   | Biochemikalien                                                                                                                                                                                                                                  |
| Product Description        | ICP-MS Tuning Solution 2 in 2% HNO <sub>3</sub> Cerium, Lithium, Thallium, Yttrium in 2% HNO <sub>3</sub> . 12 months expiry date. Traceable to NIST 31XX series. ISO 9001:2015 certified, ISO/IEC 17025:2017 and ISO 17034:2016 accredited.... |