

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Bovine Serum NMB-NBS

|                            |                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Bovine Serum                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Biozol Catalog Number      | NMB-NBS                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Supplier Catalog Number    | NBS                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Alternative Catalog Number | NMB-NBS                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Manufacturer               | NordicMubio                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Category                   | Zellen/Zellkultur                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Species Reactivity         | Bovine                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Product Description        | As a qualitative serum with an optimized protein content to study the normal protein composition of bovine serum. To identify or quantitate a normal serum protein component using a variety of immunodiffusion techniques including immunoelectrophoresis...                                                                                                                                                                                          |
| Buffer                     | Delipidated, lyophilized, stable reference serum. No preservative added.                                                                                                                                                                                                                                                                                                                                                                               |
| Formula                    | Delipidated, lyophilized, stable reference serum. No preservative added.                                                                                                                                                                                                                                                                                                                                                                               |
| Application Notes          | Normal serum protein reference sera are shipped lyophilized at ambient temperature. After arrival store at +4C. Prolonged storage may be at -20C. After opening the vial, the lyophilized content is reconstituted by adding 1 ml of sterile distilled water, mixed gently by inversion until complete dissolution is obtained. Allow to stand at ambient temperature for 5-10 minutes to reach equilibrium. Reconstituted serum may be stored frozen. |