

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Goat anti Swine haptoglobin, Clone: [Polyclonal], Monoclonal NMB-GASW/HP

Article Name	Goat anti Swine haptoglobin, Clone: [Polyclonal], Monoclonal
Biozol Catalog Number	NMB-GASW/HP
Supplier Catalog Number	GASw/HP
Alternative Catalog Number	NMB-GASW/HP
Manufacturer	NordicMubio
Host	Goat
Category	Antikörper
Application	IP
Species Reactivity	Porcine
Product Description	The reactivity of the antiserum is restricted to haptoglobin. In immunoelectrophoresis and radial immunodiffusion (Ouchterlony), using various antiserum concentrations against fresh normal swine plasma a single precipitin line is obtained which shows...
Clonality	Monoclonal
Clone Designation	[Polyclonal]
UniProt	Q8SPS7
Buffer	Delipidated, heat inactivated, lyophilized, stable whole antiserum No preservative added Total protein and IgG concentrations in the antiserum are comparable to those of pooled normal goat serum. No foreign proteins added. Reconstitute the lyophilized ant

Source	Haptoglobin is an acute phase protein synthesized in the liver. In immunoelectrophoresis it has an alpha-2 mobility. It forms a group of heterogeneous proteins consisting of four peptide chains, two alpha and two beta chains. The molecular weight is about
Application Notes	Precipitation assays. In immunoelectrophoresis use 2 ml or equivalent against 120 ml antiserum. In double radial immunodiffusion (Ouchterlony) use a rosette arrangement with 10 ml in a 3 mm diameter center well and 2 ml serum samples (neat and serially diluted) in 2 mm diameter peripheral wells.