

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Goat anti Swine Serum proteins, Clone: [Polyclonal], Monoclonal NMB-GASW/TSP

Article Name	Goat anti Swine Serum proteins, Clone: [Polyclonal], Monoclonal
Biozol Catalog Number	NMB-GASW/TSP
Supplier Catalog Number	GASw/TSP
Alternative Catalog Number	NMB-GASW/TSP
Manufacturer	NordicMubio
Host	Goat
Category	Antikörper
Product Description	In immunoelectrophoresis against pooled serum and enriched serum proteins fractions precipitation can be observed of not less than 12 individual proteins components. In precipitating techniques as immunoelectrophoresis and radial immunodiffusion (Ouc...
Clonality	Monoclonal
Clone Designation	[Polyclonal]
Buffer	Delipidated, heat inactivated, lyophilized, stable whole antiserum. 0.1% Sodium azide Reconstitute the lyophilized antiserum by adding 1 ml sterile distilled water.
Source	Pooled whole swine serum and partly purified serum fractions. Freund's complete adjuvant is used in the first step of the immunization procedure.
Formula	Delipidated, heat inactivated, lyophilized, stable whole antiserum. 0.1% Sodium azide

Application Notes	Precipitation assays. In immunoelectrophoresis use 2 µl serum or equivalent against 120 µl antiserum. In double radial immunodiffusion (Ouchterlony) use a rosette arrangement with 10 µl antiserum in 3 mm diameter center well and 2 µl serum samples (neat and serially diluted in 2 mm diameter peripheral wells).
-------------------	--