

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Normal Goose IgY (egg-derived), conjugated with FITC, Fluorescein, Aves NMB-GIGY-F

|                            |                                                                                                                 |
|----------------------------|-----------------------------------------------------------------------------------------------------------------|
| Article Name               | Normal Goose IgY (egg-derived), conjugated with FITC, Fluorescein, Aves                                         |
| Biozol Catalog Number      | NMB-GIGY-F                                                                                                      |
| Supplier Catalog Number    | GlgY-F                                                                                                          |
| Alternative Catalog Number | NMB-GIGY-F                                                                                                      |
| Manufacturer               | NordicMubio                                                                                                     |
| Host                       | Aves                                                                                                            |
| Category                   | Antikörper                                                                                                      |
| Application                | IHC                                                                                                             |
| Conjugation                | Fluorescein                                                                                                     |
| Product Description        | Normal Goose IgY (egg-derived), conjugated with FITC, Aves 0,5 mg...                                            |
| Concentration              | See vial for concentration.                                                                                     |
| Buffer                     | Phosphate Buffered Saline, 50% Glycerol with 0.075% sodium azide                                                |
| Source                     | IgY was purified from the egg yolks of unimmunized geese. The IgY was conjugated to fluorescein isothiocyanate. |
| Purity                     | Egg yolk-derived purified IgY (non affinity purified)                                                           |
| Form                       | Fluorescein conjugated                                                                                          |
| Formula                    | Phosphate Buffered Saline, 50% Glycerol with 0.075% sodium azide                                                |

#### Application Notes

Optimal working dilutions should be determined for your particular assay conditions.