

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Pooled mouse ascites without immunoglobulins NMB-NOR/AA

Article Name	Pooled mouse ascites without immunoglobulins
Biozol Catalog Number	NMB-NOR/AA
Supplier Catalog Number	NOR/AA
Alternative Catalog Number	NMB-NOR/AA
Manufacturer	NordicMubio
Category	Zellen/Zellkultur
Product Description	Ascites obtained from mice bearing different hybridomas, pooled and deleted of all immunoglobulins. It is used as a negative control for ELISA and for indirect immunofluorescence and immunoperoxidase methods, replacing the specific ascites in the tes...
Buffer	Delipidated, heat inactivated, lyophilized, stable pooled ascites. No preservative added.
Formula	Delipidated, heat inactivated, lyophilized, stable pooled ascites. No preservative added.
Application Notes	The lyophilized ascites is shipped at ambient temperature. After arrival store at +4C. Prolonged storage may be at -20C. After opening the vial, the lyophilized content is reconstituted by adding 1 ml of sterile distilled water, mixed gently by inversion until complete dissolution is obtained. Allow to stand at ambient temperature for 5-10 minutes to reach equilibrium. Reconstituted ascites may be stored frozen.